Sådan bruges det aktive, selvvarmende tæppe BARRIER® EasyWarm®/BARRIER® EasyWarm®+

Riv pakken op ved det forskårne punkt i det øverste venstre hjørne. Tag tæppet ud af pakningen for at eksponere BARRIER® EasyWarm® / BARRIER® EasyWarm®+ for ilt.

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Mölnlycke®

Referencer:

1. Torossian A, Andrzejowski J, Raeder J. A new active self-warming blanket and forced-air warming are equally effective in preventing hypothermia in mid-duration surgery: a multinational non-inferiority trial. Poster presented at: the Anesthesiology 2014 Annual Meeting 11-15 October 2014, New Orleans, LA.
2. Torossian A, Andrzejowski J, Raeder J. A new active self-warming blanket and forced-air warming are equally effective in preventing hypothermia in mid-duration surgery: a multinational non-inferiority trial. Poster presented at: the Anesthesiology 2014 Annual Meeting 11-15 October 2014, New Orleans, LA.
3. Van de Velde M, Unenge T. An Open-Label, Multicenter, Randomised, Controlled Study to Evaluate the Efficacy and Safety of the BARRIER® EasyWarm Active Self-Warming Blanket Used for Continuous Active Warming to Prevent General Anaesthesia Induced Hypothermia During the Perioperative Surgical Period. Mölnlycke Clinical Investigation Report PD-469972 Rev 01 2013. MD12-001.
4. Torossian A, Van Gerven E, Geertsens K, Horn B, Van de Velde M, Raeder J. Active perioperative patient warming using a self-warming blanket (BARRIER EasyWarm) is superior to passive thermal insulation: a multinational, multicenter, randomized trial. *Journal of Clinical Anesthesia*. 2016;34:547-54
5. Raeder J, Geertsens K, Van de Velde M, Van Gerven E, Horn B, Torossian A. Reduced hypothermia and improved patient thermal comfort by perioperative use of a disposable active self-warming blanket. A randomized multicenter trial. MD12-001. Poster presented at: 67th Annual PostGraduate Assembly in Anesthesiology; 2013 Dec 13-17; New York, NY. 2013.
6. Brandes IF, Müller C, Perl T, Russo SG, Bauer M, Bräuer A. Efficacy of a novel warming blanket: Prospective randomized trial. *Anaesthesist*. 2013;62(2):137-42.
7. Dostálová V, Schreiberova J, Bartoš M, Česák T, Habalová J, Dostálová V, et al. Thermal management in patients undergoing elective spinal surgery in prone position - A prospective randomized trial. *Ceska a Slovenska Neurologie a Neurochirurgie*. 2017;80(5):553-60.
8. Thapa HP, Kerton AJ, Peyton PJ. Comparison of the EasyWarm® self-heating blanket with the Cocoon forced-air warming blanket in preventing intraoperative hypothermia. *Anaesthesia and Intensive Care*. 2019;47(2):169-74.
9. Rosenkilde C, Vamosi M, Lauridsen JT, Hasfeldt D. Efficacy of Prewarming With a Self-Warming Blanket for the Prevention of Unintended Perioperative Hypothermia in Patients Undergoing Hip or Knee Arthroplasty. *Journal of perianesthesia nursing : official journal of the American Society of PeriAnesthesia Nurses*. 2017;32(5):419-28.
10. Koc BB, Schotanus MGM, Kollenburg JPAPAC, Janssen MJA, Tijssen F, Jansen EJP. Effectiveness of Early Warming With Self-Warming Blankets on Postoperative Hypothermia in Total Hip and Knee Arthroplasty. *Orthopedic nursing*. 2017;36(5):356-60.

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Mölnlycke, BARRIER og EasyWarm varemærker, navne og logoer er registreret globalt til en eller flere af virksomhederne i Mölnlycke Health Care gruppen. ©2022 Mölnlycke Health Care AB.